

CULTURE MEDIA FOR COSMETICS INDUSTRY ACCORDING TO ISO STANDARD

PN-EN ISO 11930:2012 Evaluation of the antimicrobial protection of a cosmetic product

NEUTRALIZERS

3102	Eugon LT 100 Broth	bottle	500,200,100,90 ml
6157	Eugon LT 100 Broth	tube	10,9 ml

CULTURE MEDIA

Identification of *A.brasiiliensis*

9007	Potato Dextrose Agar	powder	500 g
9007	Potato Dextrose Agar	bottle	500,200,100 ml
9007	Potato Dextrose Agar	plate	90 mm

Identification of *C.albicans*

9010	Sabouraud Dextrose Agar	powder	500 g
9010	Sabouraud Dextrose Agar	bottle	500,200,100 ml
9010P	Sabouraud Dextrose Agar	plate	90 mm

Colony Count

9012	Tryptic Soy Agar	powder	500 g
9012	Tryptic Soy Agar	bottle	500,200,100 ml
9012	Tryptic Soy Agar	plate	90 mm

Reference strains Epower™

0392E6	<i>Aspergillus brasiliensis</i> ATCC 16404	average number of : 1.0-9.9E+06 CFU in pellet
0443E7	<i>Candida albicans</i> ATCC 10231	average number of: 1.0-9.9E+07 CFU in pellet
0483E7	<i>Escherichia coli</i> ATCC 8739	average number of: 1.0-9.9E+07 CFU in pellet
0484E7	<i>Pseudomonas aeruginosa</i> ATCC 9027	average number of: 1.0-9.9E+07 CFU in pellet
0485E7	<i>Staphylococcus aureus</i> ATCC 6538	average number of: 1.0-9.9E+07 CFU in pellet

Detection of *Staphylococcus aureus* according to PN-EN ISO 22718

175	Maximum Recovery Dilluent	powder	500g
3531	Maximum Recovery Dilluent	bottle	500,200,100 ml
6315	Maximum Recovery Dilluent	tube	9 ml
9012	Tryptic Soy Agar	powder	500 g
9012	Tryptic Soy Agar	bottle	500,200,100 ml
9012	Tryptic Soy Agar	plate	90 mm
3102	Eugon LT 100 Broth	bottle	500,200,100,90 ml
6157	Eugon LT 100 Broth	tube	9 ml
6157	Eugon LT 100 Broth	tube szklana	10 ml
131	Baird Parker Agar base	powder	500 g
X085	Egg emulsion with pottasium tellurite	bottle	100 ml
3082	Baird Parker Agar base	bottle	500,200,100 ml
1320	Baird Parker Agar	plate	90 mm
M43	Microgen Staph Quick detection of <i>Staphylococcus aureus</i> .	100 test	lateks
M433	Microgen Staph Quick detection of <i>Staphylococcus aureus</i> .	500 test	lateks
60	Rabbit serum	vial	10x2 ml
61	Rabbit serum	vial	5 ml
3155	Catalase test	bottle	30 ml

Detection of *Pseudomonas aeruginosa* according to PN-EN ISO 22717

175	Maximum Recovery Dilluent	powder	500g
3531	Maximum Recovery Dilluent	bottle	500,200,100 ml
6315	Maximum Recovery Dilluent	tube	9 ml
9012	Tryptic Soy Agar	powder	500 g
9012	Tryptic Soy Agar	bottle	500,200,100 ml

9012	Tryptic Soy Agar	plate	90 mm
3102	Eugon LT 100 Broth	bottle	500,200,100,90 ml
6157	Eugon LT 100 Broth	tube	9 ml
6187	Eugon LT 100 Broth	tube szklana	9 ml
9001	Cetrymide Agar	powder	500 g
9001BT500	Cetrymide Agar	bottle	500,200,100 ml
9001PD90	Cetrymide Agar	plate	90 mm
MID 61 g 181804	Oxidase Strips Oxidase Strips	strips strips	50pcs 25pcs
7007B	Pseudomonas Agar P	bottle	500,200,100 ml
7007	Pseudomonas Agar P	plate	90 mm

Detection of *Escherichia coli* according to PN-EN ISO 21150

175	Maximum Recovery Dilluent	powder	500g
3531	Maximum Recovery Dilluent	bottle	500,200,100 ml
6315	Maximum Recovery Dilluent	tube	9 ml
9012D	Tryptic Soy Agar	powder	500 g
9012	Tryptic Soy Agar	bottle	500,200,100 ml
9012	Tryptic Soy Agar	plate	90 mm
3102	Eugon LT 100 Broth	bottle	500,200,100,90 ml
6157	Eugon LT 100 Broth	tube	9 ml
6187	Eugon LT 100 Broth	tube szklana	9 ml
9003	MacConkey Agar	powder	500 g
9003	MacConkey Agar	bottle	500,200,100 ml
9003	MacConkey Agar	plate	90 mm
104	Levine Eosin Methylene Blue Agar (LEMB)	powder	500 g
3033	Levine Eosin Methylene Blue Agar (LEMB)	bottle	500,200,100 ml
1011	Levine Eosin Methylene Blue Agar (LEMB)	plate	90 mm

Detection of *Candida albicans* according to PN-EN ISO 18416

175	Maximum Recovery Dilluent	powder	500g
3531	Maximum Recovery Dilluent	bottle	500,200,100 ml
6315	Maximum Recovery Dilluent	tube	9 ml
3102	Eugon LT 100 Broth	bottle	500,200,100,90 ml
6157	Eugon LT 100 Broth	tube	9 ml
6187	Eugon LT 100 Broth	tube szklana	9 ml
9013	Tryptic Soy Broth	powder	500 g
9013	Tryptic Soy Broth	bottle	500,200,100 ml
261	Letheen Broth base	powder	500 g
3225	TWEEN 80	bottle	200,100ml
3541	Letheen Broth z TWEEN 80	bottle	500,200,100 ml
9010	Sabouraud Dextrose Agar	powder	500 g
9010	Sabouraud Dextrose Agar	bottle	500,200,100 ml
9010	Sabouraud Dextrose Agar	plate	90 mm
9007	Potato Dextrose Agar	powder	500 g
9007	Potato Dextrose Agar	bottle	500,200,100 ml
9007	Potato Dextrose Agar	plate	90 mm
114	Sabouraud Dextrose Agar z chloramfenikolem	powder	500 g
3020	Sabouraud Dextrose Agar z chloramfenikolem	bottle	500,200,100 ml
1231	Sabouraud Dextrose Agar z chloramfenikolem	plate	90 mm
7008	Corn meal Agar Base	powder	500 g
3225	TWEEN 80	bottle	200ml
3225	TWEEN 80	bottle	100ml
7008	Corn meal Agar z 1% Tween80	bottle	500,200,100 ml
7008	Corn meal Agar z 1% Tween80	tube	90mm
7000	Potato Dextrose Agar z chloramfenikolem	bottle	500,200,100 ml
7000	Potato Dextrose Agar z chloramfenikolem	plate	90 mm
998800QCA	Gram PVP Kit	bottle	4x250ml

